

Brown Stacking Commode

UTO1021-2, UTO1022-2



USER MANUAL



WARNING: Please read this Manual carefully before using the **Brown stacking commode (with and without arms)** as it contains important information regarding the safe use of the product.

Contents

1.	Introduction	4
2.	Description of product	4
3.	Instructions for use	4
3.1	Adjust commode height	4
3.2	Detach commode arms – UTOC1022-2 only	5
3.3	Removal of the commode buckets	5
4.	Cleaning	6
5.	Care and maintenance	6
6.	Technical Specification	7
7.	Warranty	8
8.	Disposal advice	9
9.	Life expectancy.....	9
10.	Contact Us	10

1. Introduction

These instructions aid the user and/or carer about the safe and efficient use of this product, including assembly, cleaning and maintenance, and operation of the product.

2. Description of product

The Ultimate Healthcare commodes are designed to assist those who have difficulty in accessing or using a standard toilet. It can also be used as a temporary chair. All Ultimate Healthcare commodes have a hinged plastic toilet seat, commode bucket with lid and vinyl lift off seat.

3. Instructions for use

Remove vinyl seat to use commode. Replace vinyl seat to use commode as a temporary chair.

3.1 Adjust commode height

- To adjust the height, remove each 'E' clip in turn. Slide each leg to the required height.
- Align the holes and replace 'E' clip.
- Ensure that all legs are secured at the same position to prevent rocking.



WARNING:

- Do not extend legs in a manner so that the commode is slanting. This will compromise the safety of the user.

3.2 Detach commode arms – UTOC1022-2 Only

- To detach arms for transfer, depress the metal buttons each side of the commode arms fixing sockets.
- Using both hands, remove commode arm by lifting arm from fixing socket.
- To refit arms, locate arms into fixing sockets and push down to stop.
- Depress buttons on each end of arm, ensuring buttons “click” firmly into place.

3.3 Removal of the commode buckets

To remove the commode bucket, simply lift the seat of the toilet and fix the lid provided onto the bucket to avoid splashing and carefully raise the bucket and dispose of the contents into a proper sluice or WC.



WARNING:

- Do not use an open drain, sink or basin. This could cause a serious health risk.

Clean out the residue with a suitable disinfectant and refit to the commode after use.

4. Cleaning

To clean vinyl seats use a suitable detergent and warm water, wipe seats and dry with disposable paper towel. Clean steel frame and toilet seat with a suitable detergent and warm water, wash framework and arm rests thoroughly paying particular attention to the toilet seat hinges and the underside of the toilet seat.

Can be autoclaved at 80°C.

Allow equipment to dry prior to stacking and storage.



CAUTION:

- Under no circumstances should the commode be moved whilst occupied. It is advised that you check the frame of the product annually for signs of wear and tear, dimpled or flaking paint or cracks in the frame should be reported to your supplier immediately. Should there be any such evidence, please refrain from using the commode.

5. Care and Maintenance

This product is maintenance free so long as the cleaning and safety instructions shown above are adhered to.

Check the product regularly for damage and to ensure that it is securely assembled.

6. Technical specification

Frame is composed of welded, Powder coated mild steel tube.

The seat pad is high frequency welded PVC with a foam and wood filling.

The extension legs are made from chrome plated mild steel tube.

	Brown Stacking Commode	Brown Stacking Commode with removeable arms
Product code	UTOC1021-2	UTOC1022-2
Seat pad width	440mm	
Seat pad depth	440mm	440mm
Seat height (inc seat pad)	490mm – 640mm	
Width between armrests	505mm	
Max. height	930mm	
Min. height	780mm	
Max. width	550mm	
Min. width	550mm	
Max. depth	550mm	
Min. depth	520mm	460mm
Max user weight	170kg/ 26st	
Weight of product	9kg	

7. Warranty

The Brown Stacking Commode Chairs are warranted for 12 months from the date of purchase against manufacturing or design defects.

The warranty is subject to the following conditions:

- It does not apply to wear and tear, inappropriate or incorrect use.
- Defective set up by the user or third parties. Negligent or careless treatment.
- Use of unsuitable or unapproved spares or unauthorised alteration.
- The warranty does not extend to consequential costs resulting from fault clearance, loss of earnings or expenses

This in no way affects your statutory rights.

8. Disposal advice

This product must be thoroughly cleaned/sterilised prior to disposal. The product and /or its components are of standard non-hazardous materials therefore can be disposed of without any potential environmental hazards providing that local recycling arrangements are sought and adhered to.

9. Life expectancy

We estimate a life expectancy of 2 years for this product, provided that:

- It is used in strict accordance with the intended use as set out in this document.
- All service and maintenance requirements are met.

The estimated life expectancy can be exceeded if the product is carefully used and properly maintained.

The life expectancy can also be considerably reduced by extreme or incorrect usage.

The fact that we estimate a life expectancy for this product does not constitute an additional warranty.

10. Contact Us

Ultimate Healthcare customer service is available to answer questions from 8.30am to 5.00pm, Monday to Friday for business and upgrades. Our sales team is available for pricing and order processing and status queries. Our Technical Support team is available to answer your queries related to troubleshooting, please contact the following:

Ultimate Healthcare Limited

Nutsey Lane
Calmore Industrial Estate
Totton
Southampton
SO40 3XJ

Tel: 0333 321 8996

Email: Enquiries@ultimatehealthcare.co.uk

Please Note: Information in this user manual is subject to change without notice and does not represent commitment on the part of Ultimate Healthcare Limited. No part of this manual may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying and recording for any purpose without the written permission of Ultimate Healthcare.

Page left intentionally blank

Ultimate Healthcare Limited.
Calmore Industrial Estate
Totton, Southampton
SO40 3XJ
United Kingdom



Doc no.
UTOC1021-2.1022-2.IFU-3

Issue 3.0
May 2025

www.ultimatehealthcare.co.uk