

Overbed Table

UHOT1010-2 – No Castors

UHOT1020-2 – With Castors

USER MANUAL





WARNING: Please read this Manual carefully before using the overbed table as it contains important information regarding the safe use of the product.

Contents

1.	Description of product	4
2.	Assembly Instructions	4
3.	Cleaning Instruction	6
4.	Technical Specifications	6
5.	Warranty	6
6.	Contact Us	7

1. Description of product

This over bed table has a laminated ash top and two independent points for adjusting the height, which allows the tabletop to be angled. Edging strips on both sides of the tabletop ensure items don't roll off when it is used in the angled position, which can be an ideal position for ensuring a comfortable read, crossword or hobby. Laid flat the table tops good sized surface is perfect for meals or snacks. The four castors (**Product UHOT1020-2**) allow the table to easily glide to and from the bed or chair, ideal for those with limited strength or mobility.

2. Assembly Instructions

- 1) Attach castors by locating pin in the hole at the bottom of the table base frame and pushing them in to place, do not hit (**Product UHOT1020-2 only**)
- 2) Turn base tabletop over and insert "H" frame into the holes at the top with label facing up
- 3) Secure H frame using "E" clips
- 4) Insert silver L shape pieces into the back of the tabletop with spring pins fully inserted through both black brackets
- 5) Flip tabletop over and insert silver L shape pieces into the "H" frame
- 6) Secure using black hand screws to the desired height and angle



WARNING:

The table is not load bearing and should not be used to support objects over 10kg (22 lbs). Please do not allow people to lean on the table.



CAUTION:

Periodically check to ensure all components are present, correctly assembled and suitably tightened.

NOTE:

To adjust the height and/or angle of the table, first loosen the 2 shorthand screws by turning anticlockwise. Adjust the table to suit, and secure by turning hand screws clockwise.

3. Cleaning Instruction

Use a non-abrasive mild detergent along with warm water & a soft cloth to clean the table. Rinse with clean warm water & a soft dry cloth.



CAUTION:

- It is not recommended to use any bleach, very hot water or any abrasive surface cleaners, as this may result in damage or mark the table surface.

4. Technical specification

	UHOT1020-2	UHOT1010-2
Maximum load weight:	10kg/ 22 lbs	10Kg/22lb
Length:	697 mm	680mm
Width:	413 mm	413mm
Height (From floor to top of table):	680 – 980 mm	630 – 930 mm

5. Warranty

This product is guaranteed to be free of defects in materials or workmanship for one year from date of purchase.

This product has been built to high standards of workmanship and has been carefully inspected prior to shipment. This warranty is an expression of the confidence we have in the workmanship and materials used in our products. In the unlikely event of there being a defect found with the product, please contact your supplier and they will either replace the product or supply parts for the product to be repaired at their discretion.

This warranty does not cover failure of the device caused by misuse, negligence or fair wear and tear of the product.

6. Contact Us

Ultimate Healthcare customer service is available to answer questions from 8.30am to 5.00pm, Monday to Friday for business and upgrades. Our sales team is available for pricing and order processing and status queries. Our Technical Support team is available to answer your queries related to troubleshooting, please contact the following:

Ultimate Healthcare Limited

Nutsey Lane
Calmore Industrial Park
Totton
Southampton
SO40 3XJ

Tel: 0333 321 8996

Email: enquiries@ultimatehealthcare.co.uk



Please Note: Information in this user manual is subject to change without notice and does not represent commitment on the part of Ultimate Healthcare Limited. No part of this manual may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying and recording for any purpose without the written permission of Ultimate Healthcare

Ultimate Healthcare Limited.
Calmore Industrial Estate
Totton, Southampton
SO40 3XJ
United Kingdom



Doc no.
UHOT10..-2.IFU

Issue 3.0
September 2025