

## Carefree Overlay

Pressure reducing foam mattress overlay



The **Carefree Overlay** is used on top of a standard mattress. It comprises of a castellated foam to help provide pressure reduction.

The mattress overlay is ideal for patients considered to be up to a **medium risk** of pressure ulcer development and can be used in a **community and long term care setting**.

Available in both 2" and 4" height versions as a **single** or **double** size to provide flexible care options.

### Features & Benefits



#### Castellated foam

Offering excellent pressure reduction, comfort and support.



#### Non-turn Overlay mattress

Not requiring turning or rotating significantly minimising manual handling risks.

#### 2-way stretch, waterproof, vapour permeable PU cover

Improving patient comfort and maintaining skin integrity by reducing shear and friction forces whilst aiding infection control.

For more information on this product please contact us on tel: 0333 321 8996 or email: [sales@ultimatehealthcare.co.uk](mailto:sales@ultimatehealthcare.co.uk)

[www.ultimatehealthcare.co.uk](http://www.ultimatehealthcare.co.uk)



## Specification

<b>Product code:</b>	UPRF3478-2	<b>Dimensions:</b>	86 x 198 x 5cm	(34 x 78 x 2") <b>Single</b>
	UPRF3478		86 x 198 x 10cm	(34 x 78 x 4") <b>Single</b>
	UPRF5478-2		137 x 198 x 5cm	(54 x 78 x 2") <b>Double</b>
	UPRF5478		137 x 198 x 10cm	(54 x 78 x 4") <b>Double</b>

*Additional sizes available on request*

<b>Max user weight:</b>	135 kg / 21 stone (for 2" Overlay)	<b>Risk category:</b>	Up to medium
	165 kg / 26 stone (for 4" Overlay)	<b>Product weight:</b>	3kg - 8.5kg
		<b>Fire retardancy:</b>	BS7177 (Crib 5 Composite)
		<b>Warranty:</b>	2 years


**Cleaning:** Can withstand up to 10,000 ppm chlorine (to be used in conjunction with local infection control guidelines) (can be laundered up to 84°C)

## Pressure Mapped Patient

**Age:** 44      **Weight:** 83kg  
**Gender:** Male      **Height:** 170cm



2" Overlay



4" Overlay

## Accessories

Carefree overlay replacement covers

<b>UPRF3478-2C</b>	86 x 198 x 5cm	for UPRF3478-2
<b>UPRF3478C</b>	86 x 198 x 10cm	for UPRF3478
<b>UPRF5478-2C</b>	137 x 198 x 5cm	for UPRF5478-2
<b>UPRF5478C</b>	137 x 198 x 10cm	for UPRF5478

## Options available

Carefree replacement mattresses

<b>UPRF3678</b>	91 x 198 x 15cm	Single
<b>UPRF5478-6</b>	137 x 198 x 15cm	Double

For more information on this product please contact us on tel: 0333 321 8996 or email: [sales@ultimatehealthcare.co.uk](mailto:sales@ultimatehealthcare.co.uk)

[www.ultimatehealthcare.co.uk](http://www.ultimatehealthcare.co.uk)

