

Tamora II

Dynamic Mattress Overlay



The Tamora II is a lightweight, portable mattress overlay system that is suitable for use in a wide range of care settings from general ward and acute care to private or nursing home care. Ideal for patients considered to be up to a high risk of pressure damage the Tamora II is an optimal pressure care solution.



Lightweight, portable and easy to use

The Tamora II is a lightweight, portable system which is designed to be used on a patient's existing mattress or divan bed and can be secured safely in place by **elastic securing straps** which are affixed to the base of the overlay.

Simple and concise control panel operations allow optimised patient care with a **Comfort Control dial** which allows nursing staff to adjust pressures to suit individual patient requirements and provide enhanced comfort and support whilst maintaining a good level of protection and therapy.



Assistance with infection control procedures

The Tamora II has been designed to ease the burden on nursing staff in the implementation of cleaning and decontamination procedures.

Key features such as an easy to clean, waterproof cover with **welded seams**, **enclosed zipper** and integral **covered air pipe** will assist hospital staff with infection control initiatives and minimise the risk of cross infection.



Safety and peace of mind

Patient and carer safety are key factors in the system design with many features which specifically enhance the safety of patients and caregivers alike.

An easy to use **cable management system** located on both sides of the mattress cover eliminates tripping hazards from loose cabling and reduces potential for power damage cable. While a timed **control panel lockout** feature prevents unintentional activation of any functions which affect the patient, after a period of inactivity the control panel will lockout.

Features

Power Unit

- Comfort control dial
- 7 litres per minute air flow
- Alternating 2-cell cycle
- Variable cycle times (10/15/20 mins)
- Automatic control panel lockout
- Timed static mode
- Transport mode - up to 24 hours
- Audible and visual alarms
- Adjustable hanging brackets


Mattress

- Static head cells
- CPR Valve
- Strong knitted cell retention bands
- Maximum patient weight 200kg / 31.5 stone

Cover

- Soft, waterproof vapour permeable PU cover
- 360 degree zipper
- Welded seams and enclosed zip
- Integral covered air pipe
- Cable management system
- Elastic securing straps

Tamora II Overlay Mattress - Pressure Mapping Patient: Male at a height of 5ft 7, weighing 82kg



Technical Specification

Product code:	UPRA3478-2	Cycle time:	10 /15 / 20 mins (2-cell cycle)
Risk category:	Up to a high risk	Pressure range:	30-70mmHg
Max user weight:	200 kg / 31.5 stone	Supply voltage:	230V, 50 Hz
Mattress dimensions:	2000mm x 850mm x 130mm	Fuse rating:	T1A
Power unit dimensions:	320mm x 130mm x 230mm	Current rating:	0.12 A
Product weight:	4.9kg Mattress / 3.5kg Power Unit	Classification:	Class II, Type BF
Cell material:	Nylon PU	Compliance:	EN 60601-1, EN 60601-1-2
Cover material:	Polyester PU	Warranty:	2 years

Cleaning: Can withstand up to 10,000 ppm chlorine (to be used in conjunction with local infection control guidelines) (can be laundered up to 84°C)

Ultimate Healthcare,

Calmore Industrial Estate, Nutsey Lane, Totton, Hampshire, SO40 3WW

tel: 0333 321 8996, fax: 0845 051 0035, email: sales@ultimatehealthcare.co.uk

www.ultimatehealthcare.co.uk