


## Trinity Plus Premier

Visco elastic foam mattress replacement



The **Trinity Plus Premier** replaces a standard mattress. It comprises of a castellated visco elastic foam insert to provide pressure redistribution by conforming to a users weight and temperature. The mattress also features a reinforced base and side wall.

The mattress replacement is ideal for patients considered to be up to **very high risk** of pressure ulcer development and can be used in a **community and long term care** setting.


Available in **single** and **double** sizes to provide flexible care options.

### Features & Benefits



#### Castellated visco elastic foam

Allows for immersion, comfort and postural support whilst re-distributing pressure away from high risk areas.



#### Reinforced side walls and base

Offering stability and support whilst assisting patients and their caregivers with transfers.



#### Non-turn mattress replacement

Not requiring turning or rotating significantly minimises manual handling risks.

#### 2-way stretch, waterproof, vapour permeable PU cover

Improving patient comfort and maintaining skin integrity by reducing shear and friction forces whilst aiding infection control.

High frequency welded seams and a fully concealed zip which reduces the possibility of fluid ingress to aid infection control.

For more information on this product please contact us on tel: 0333 321 8996 or email: [sales@ultimatehealthcare.co.uk](mailto:sales@ultimatehealthcare.co.uk)

[www.ultimatehealthcare.co.uk](http://www.ultimatehealthcare.co.uk)






## Specification

<b>Product code:</b>	UPRVP3678	<b>Dimensions:</b>	91 x 198 x 16cm	(36 x 78 x 6") <b>Single</b>
	UPRVP5478-6		137 x 198 x 16cm	(54 x 78 x 6") <b>Double</b>
<b>Max user weight:</b>	260 kg / 41 stone		<i>Additional sizes available on request</i>	
<b>Risk category:</b>	Up to very high	<b>Fire retardancy:</b>	BS7177 (Crib 5 Composite)	
<b>Product weight:</b>	13kg	<b>Warranty:</b>	2 years	

**Cleaning:** Can withstand up to 10,000 ppm chlorine (to be used in conjunction with local infection control guidelines) (can be laundered up to 84°C)

## Pressure Mapped Patient

<b>Age:</b> 44	<b>Weight:</b> 83kg
<b>Gender:</b> Male	<b>Height:</b> 170cm



6" Mattress replacement

## Accessories

### Mattress extension squab infill

<b>USPF3408</b>	for use with	UPRVP3678
<b>USPF3406-KIT</b>	for use with	UPRVP3678 on a Cura Community bed

### Trinity Plus Premier mattress replacement covers

<b>UPRVP3678C</b>	91 x 198 x 16cm	for UPRVP3678
<b>UPRVP5478-6C</b>	137 x 198 x 16cm	for UPRVP5478-6



For more information on this product please contact us on tel: 0333 321 8996 or email: [sales@ultimatehealthcare.co.uk](mailto:sales@ultimatehealthcare.co.uk)

[www.ultimatehealthcare.co.uk](http://www.ultimatehealthcare.co.uk)